

DDL Test Record 10-10-2007

| | |
|---|--------------------|
| Test Details | |
| Date of test | 2007-10-10 |
| Tester Name | Caroline D. Cooper |
| Version of DDL | DDL version 2091M |
| Build date | Oct 1, 2007 |
| Browser Type (Firefox, IE) and version | Firefox 2.0.0.7 |
| OS type and version (Windows XP, Vista) | Windows XP |

| Test | P/F/NR | Details (include Mantis bug number) |
|--|--------|---|
| Uninstall | Pass | Followed instructions, did not run into problems |
| Install | Pass | settings link is broken, used http://ci.myhealthspace.com:9080/DDL/ directly |
| Configure DDL | Fail | Update DICOM Configurations" leads to 404 error, though settings are correctly saved |
| Test 1: setting active worklist | Pass | |
| Test 2: uploading data to MedCommons via DDL | Pass | Plain upload- Pass Cancel and delete- pass Cancel and restart- pass |
| Test 3: Downloading data from MedCommons | Fail | Export to local file system and export to DICOM device both work sometimes but more often fail. Sometimes a total file size and percent completed are not displayed; in this case the file will continue to indicate it is being "downloaded" indefinitely. This is the most common problem. The other problem is that errors in configuring a remote DICOM device or a directory on the local filesystem are not detected until the download is complete, and once the download is complete and export failed there is no way to attempt export again. |
| Test 4: round trip | NR | Test is not coherent enough to run. |
| | | |