

DDL Test Record 10-12-2007

Test Details	
Date of test	2007-10-12
Tester Name	Caroline D. Cooper
Version of DDL	DDL version 2091M
Build date	Oct 1, 2007
Browser Type (Firefox, IE) and version	Firefox 2.0.0.7
OS type and version (Windows XP, Vista)	Windows XP

Test	P/F/ NR	Details (include Mantis bug number)
Uninstall	Pass	Followed instructions, did not run into problems
Install	Pass	settings link is broken, used http://ci.myhealthspace.com:9080/DDL/ directly
Configure DDL	Fail	Update DICOM Configurations" leads to 404 error, though settings are correctly saved
Launch DDL	Pass	Gives confusing messages on startup ("user must be logged in" when user is logged in to a medcommons appliance already, etc.) but resolves itself.
Test 1: setting active worklist	Pass	<ul style="list-style-type: none">- balloon pops up correctly and corresponds with About page- Uploads to the worklist indicated- AS FAR AS I CAN TELL, it uploads to the worlisk indicated when the user selects the folder to be uploaded, even if the worklist is changed before the folder finishes being imported into the DDL. I noticed this while attempting to verify "uploads to worklist indicated" as I dealt with the problems i noted below. Is this worth testing specifically?
Test 2: uploading data to MedCommons via DDL	Fail	<p>Long delay before anything shows up on Status page, makes it difficult to determine whether upload has started or not. "Clear all transactions" cleared all transactions that appeared on the status page, but did not clear those that had not yet shown up (which sometimes took over a minute to start).</p> <p>Several uploads failed, which may be because I had attempted to upload the same folder (/DICOM/Adrian on the shared drive) three times in quick succession. Uploaded /DICOM/Two Images with no problems.</p>

		<p>Upon a successful upload, it attempts to open the CCR (? it opens a new window/tab at any rate). This leads to an error page. Closing this page and opening the CCR again causes the CCR to display properly. The same DICOM can be uploaded many times. Each turns into a separate entry on the worklist.</p>
Test 3: Downloading data from MedCommons	Fail	<p>The folder "DICOM/Two image" also known as "alpha test" was uploaded twice in a row. Subsequently attempting to download the DICOM images from the resulting CCR gave the follow error message in a balloon: "Matching identifier in database MGH 1713806380 should have a single matching entry, not 2 Files not sent to DICOM CSTORE target"</p> <p>With a different CCR: Download to Osirix on different machine: Failed; Download does not start. Cancelling and restarting the download causes it to download correctly.</p> <p>Download to file: Failed; download does not start. Cancelling and restarting the download causes it to download correctly.</p>
Test 4: round trip	NR	Test is not coherent enough to run.